



R&D SOP098 Pipetting

Document Title: Pipetting

Document Number: R&D SOP098

Staff involved in development: <i>Job titles only</i>	Senior R&D Manager, R&D Operational Manager, Clinical Project Managers and Research Pathology
Document author/owner:	Senior R&D Manager
Directorate:	Research and Development
Department:	Research and Development
For use by:	NHS Staff Trust-Wide
Review due:	DATE
<p><u>THIS IS A CONTROLLED DOCUMENT</u></p> <p>Whilst this document may be printed, the electronic version maintained on the Trust's Intranet is the controlled copy. Any printed copies of this document are not controlled. ©Royal Papworth Hospital NHS Foundation Trust. Not to be reproduced without written permission.</p>	

Summary of Amendments

Version Number	Modification:

Key Points of this Document

1 Purpose and Contents

- a. To ensure accuracy and precision pipettes must be handled with care and maintained regularly.
- b. Air displacement pipetting is highly accurate for standard pipetting applications. However, conditions, such as atmospheric pressure as well as the specific gravity and viscosity of the solution, may affect the performance of air displacement pipettes.

2 Roles & Responsibilities

- a. A pipette use should only be performed by individuals fully trained by a competent member of the R and D team.
- b. Following a period of supervision (depending on the individual needs of the trainee) there will be an informal assessment.

3 Policy

- a. This SOP is mandatory and, as per the Trust's Information Governance and Records Management framework, non-compliance with it may result in disciplinary procedures.

4 Procedure

- a. Types of Pipettes
 1. Gilson and similar models
Most common pipette used in most applications.

R&D SOP098 Pipetting

Figure 1:



2. Disposable/Transfer Pipette- Pasteur Pipettes

This is the most basic type of pipette; and can be used for rough measurements only. Always use a new pipette and dispose of it after use either in a biohazard bin or a sharps bin.

- Aspirate liquid at a 90-degree angle.
- Dispense at a 45-degree angle.
- Touch off to make sure all the liquid is dispensed.

3. Finn pipette C1/pipette boy and similar model

Serological pipettes are used in cell and tissue culture applications and in general laboratory liquid dosage, when over one ml volumes are pipetted. Serological pipettes are made of glass or polystyrene. Plastic, disposable pipettes are useful in applications where sterility is

R&D SOP098 Pipetting

a requirement. The speed of both aspiration and dispensing can be adjusted separately to work with a variety of liquids. A pump or a pipette controller is used to draw up liquid. Measure solution from the bottom of the meniscus, the crescent-shaped surface of the liquid that is visible in the pipette.

Figure 2




b. Pipette Tips

1. There are both standard and multi-purpose tips for the use of laboratory applications with varying performance requirements ranging from very high accuracy to reagent dispensing with greater tolerance. Sterile standard tips are available for applications demanding the highest level of purity from pathological microorganisms.
2. Tips may also have filters which are used when processing samples for molecular work.
3. Make sure to use the correct size tip for the pipette you are using. Tips come in different sizes to accommodate the difference ranges of volumes the pipette is used for.
4. To fit a disposable pipette tip, hold the micropipette in one hand, lower the pipette into the tip, and use a slight twisting movement to seat the tip firmly on the tip holder of the micropipette to ensure an air-tight seal.
5. To manually fit a disposable tip, hold the pipette in one hand and use a slight twisting movement to seat the tip firmly on the shaft of the micropipette and to ensure an air-tight seal.
6. To avoid touching contaminated tips, hold the pipette over the waste container and press the tip ejector push-button. Eject the used tip and store the pipette in an upright position.

R&D SOP098 Pipetting

7. To protect your pipette, avoid tapping the tip onto the pipette like a hammer. Tips are available in TIPACK racks for easy mounting
 8. Always eject used tips. The residual liquid can damage your pipette.
- c. Adjusting the volume for a Gilson pipette
1. Check on top the pipette what is the range before adjusting the thumbwheel. The volume is set by turning the thumbwheel or the push-button. The pushbutton makes it easier and quicker to set volumes, especially when wearing gloves. The thumb wheel may be turned to slowly reach the required setting.
 2. To avoid internal damage to your pipette, never attempt to force the volume setting beyond the limits. If the dial is difficult to move you may have reached the limit. Check before twisting further. If a number is displayed in red the pipette is out of range.
- d. Pipetting technique- Forward pipetting

Figure 4



1. Fit the tip onto the pipette tip cone.
 2. Press the operating button to the first stop.
 3. Place the tip just under the surface of the liquid (2-3mm) and smoothly release the operating button up to the starting position.
 4. Wait one second.
 5. Carefully withdraw the tip from the liquid, touching against the edge of the container to remove excess from the outside of the tip.
- e. Cleaning of the outer surface of the pipette
1. Remove Tip-Holder and Tip-ejector.
 2. Keep the tip-ejector button depressed and grip the top of the tip-ejector with the other hand.
 3. Gently rotate the tip-ejector counter clockwise and separate its connector from the operating rod.
 4. Pull the tip-ejector away from the body of the pipette.
 5. Wipe the tip ejector with a soft-cloth or lint-free tissue lightly coated in 70% alcohol solution.

R&D SOP098 Pipetting

6. Wipe all exposed surfaces of the body, button, operating rod, and tip holder with a soft-cloth or lint-free tissue lightly coated with the alcohol solution.
7. Remove all scuff and wear marks.
8. Remove build up in concave surfaces.
9. If the pipette is very dirty, a brush with soft plastic bristles may be used.
10. Refit the tip ejector and allow the pipette to dry.
11. Wipe outer surfaces with a soft-cloth or lint-free tissue lightly coated with the 70% alcohol solution to remove any remaining residue from cleaning.
12. Allow time for alcohol to evaporate.

f. Testing performance

If you are worried about the accuracy of the pipette, use distilled water and weigh out the total volume it dispenses. Repeat this 10 times to get an average. If the pipette is not measuring accurate volumes please inform a member of staff so the pipette can be taken out of use until recalibrated.

4.2 Servicing

- a. Pipettes in use should be serviced annually.
- b. Pipettes need to be decontaminated/ cleaned with 70% alcohol or a Clinell wipe. A DN418, and a decontamination form must be filled in to give to the engineer. Include a list of pipettes (make and model) that require servicing.
- c. All tips that are being used for the pipettes must be provided to complete calibration.
- d. Electronic pipettes must be charged prior to servicing.
- e. Pipettes will be serviced when the Cambridge University in the HLRI have their pipettes serviced. If pipettes need to be serviced or repaired at a different time contact the Pipette Doctor <http://calibrationpipetterepair.com/>.

5 Health and Safety using a Pipette**5.1 COSHH**

- a. Staff members need to read all COSHH forms related to chemicals used to clean centrifuge. These include: Clinell wipes and or 70% alcohol.



R&D SOP098 Pipetting

6 Risk Management / Liability / Monitoring & Audit

- a. The R&D SOP Committee will ensure that this SOP and any future changes to this document are adequately disseminated.
- b. The R&D Department will monitor adherence to this SOP via the routine audit and monitoring of individual clinical trials and the Trust's auditors will monitor this SOP as part of their audit of Research Governance. From time to time, the SOP may also be inspected by external regulatory agencies (e.g. Care Quality Commission, Medicines and Healthcare Regulatory Agency).
- c. In exceptional circumstances it might be necessary to deviate from this SOP for which written approval of the Senior R&D Manager should be gained before any action is taken. SOP deviations should be recorded including details of alternative procedures followed and filed in the Investigator and Sponsor Master File.
- d. The Research and Development Directorate is responsible for the ratification of this procedure.

Further Document Information

Approved by: <i>Management/Clinical Group</i>	<i>Directorate</i>	Research and Development Directorate					
Approval date: <i>(this version)</i>	Current approved version date						
Ratified by Board of Directors/ Committee of the Board of Directors:	STET						
Date:	N/A						
This document supports: <i>Standards and legislation</i>	Medicines for Human Use (Clinical Trials) Regulations 2004 and all associated amendments. UK Policy Framework for Health and Social Care Research (2018)						
Key related documents:	Trust Research Policy Trust Policy DN1 Document Control Procedures						
Equality Impact Assessment: Does this document impact on any of the following groups? If YES, state positive or negative, complete Equality Impact Assessment Form available in Disability Equality Scheme document DN192 and attach.							
Groups	Disability	Race	Gender	Age	Sexual orientation	Religious & belief	Other

R&D SOP098 Pipetting

Version 1.0 Review Date: April 2026

Page 7 of 9



R&D SOP098 Pipetting

Yes/No	NO	NO	NO	NO	NO	NO	NO
Positive/Negative							
Review date:			DATE				

I certify the contents of this SOP has been reviewed and ratified

DocuSigned by:
Patrick Calvert
B1A52758BFFF421

26-Aug-2023

.....

.....

Signed by Dr Patrick Calvert, Clinical Director of R&D

Date

SOP Release Date:



R&D SOP098 Pipetting

Appendix - Pipettes risk assessment

<p>Description of Activity Measuring accurately liquids. Often liquids are biohazardous including blood and lavage samples.</p>	<p>Significant Hazards Manual handling Biohazard- Needlestick injury</p>		
<p>Frequency: Multiple times during the day</p>	<p>Duration: variable depending on project</p>		
<p>Adverse Effects</p> <ul style="list-style-type: none"> • Repetitive strain injury (musculoskeletal injury) 	<p>People at Risk</p> <ul style="list-style-type: none"> • Research team • Clinical Trial Co-ordinator • Tissuebank staff 		
	<p>Number of People Affected: various</p>		
<p>What Precautions exist to control the risk:</p> <ul style="list-style-type: none"> • Frequency of pipette use per person is relatively low however if workload increases can purchase electronic pipettes. • Staff competency is tested to ensure pipettes are properly used. • Encourage Staff to use ejector button instead of removing pipette tips manually. • Gloves should be worn when using a pipette especially if handling biohazardous material. • All used pipette tips and Pasteur pipettes should be put in clinical waste bin or biobin when used. • In case there is skin contact with spillage or sharps injury; Occupational health (M-F 9am-5pm) or ALERT team (bleep 432/521) MUST be notified. • Complete a Decontamination certificate after cleaning a pipette that is going to be serviced. <p>Risk Rate:</p>			
<p>Likelihood 1</p>	<p>Severity 2</p>	<p>RRN 2</p>	<p>H S M L N</p>
<p>Are these arrangements satisfactory: Yes</p>			